

Anti-human phospho-231 TAU mab 7D11**Catalogue number**

847-0102009401

847-0102009403

Package size

100 µg

1 mg

For research use only.

Product:	monoclonal antibody
Clone:	7D11
Host:	mouse
Isotype:	IgG1 _k
Immunogen:	Synthetic peptide KKVAVVRT(Pi)PPKSPSS
Specificity:	human TAU phosphorylated at amino acid position 231 (T)
Purification:	purified by affinity chromatography
Purity:	> 95% by SDS-PAGE
Solution:	PBS pH 7.4 without additives
Storage:	at -20 °C (repeated thawing and freezing should be avoided)
Applications:	ELISA, Western Blot, Immunohistochemistry
Dilution:	1 µg/ml in Western Blot
Reactivity:	PHF TAU, phosphorylated TAU at amino acid position 231 (T)



Roboscreen GmbH
Hohmannstraße 7
04129 Leipzig • Germany

CEO
Dr. Ingolf Lachmann

Phone: +49 341 989 734 0
Fax: +49 341 989 734 199
info@roboscreen.com
www.roboscreen.com

Ust-IdNr. DE 813064314
County Court Leipzig
HRB 17503
Leipzig